

**MS Schippers****MS SCANGEL****5 L e****Hypo-allergenic water soluble contact medium  
for ultrasound transmission****1609269****Directions**

Apply MS Scangel directly to diagnostic area, the amount required will depend on the size of the area. It is easily removed with a towel or tissue after use.

MS Scangel is made of water-soluble polymer. It is hypo-allergenic, odorless and nonstaining. It spreads easily and evenly.

**Schippers Export EU BV All located in:****Schippers Africa BV****Schippers APAC BV****Schippers LATAM BV****Schippers ME BV**

Tel: +31(0)497-339774

Schippers Europe BV

Smaragdweg 60

5527 LB Hapert

Tel: +31(0)497-339 779

**België:**

Schippers BVBA

Grens 114

2370 Arendonk (BE)

Tel: 014-672356

**Deutschland:**

Schippers GmbH

Kölner Straße 62

47647 Kerken (D)

Tel: 02833-92360

22 LV 03

The  
**Schippers**  
Group

**France:**

Schippers France Sarl

ZA de la Retaudais - CS

10009

35137 BEDEE (FR)

Tél: 02 99 61 40 40

**España:**

Schippers Agrícola SL

Ctra. Montmeló 76

08403 Granollers

Barcelona - España

Tel. 935 689 128

**Italia:**

Schippers Italia Srl.

Via Belvedere 84

24022 Alzano

Lombardo (BG)

Tel: 035-4490369

Download app!  
Scan this product

**1609269**